

A
L
L
O
Y
G
U
I
D
E

Alloy Name	Symbol	Element	%	Colour	Use
Argelite A	AU	GOLD	2	SILVER	CAST RESTORATIONS PORCELAIN BONDING 81% precious metal
	PD	PALLADIUM	79		
	CU	COPPER	10		
	GA	GALLIUM	8.8		
Argelite 50	AU	GOLD	0.1	SILVER	CAST RESTORATIONS PORCELAIN BONDING 50% precious metal
	PD	PALLADIUM	49.75		
	AG	SILVER	40		
	IN	INDIUM	5		
	SN	TIN	4		
	ZN	ZINC	1		
	RU	RUTHENIUM	TRACE		
Argeloy NP	CO	COBALT	59.5	SILVER	CAST RESTORATIONS PORCELAIN BONDING Nickel free non-precious metal Suitable for NHS Cast Crowns Suitable for NHS Bonded Crowns
	CR	CHROMIUM	31.5		
	MO	MOLYBDENUM	5		
	SI	SILICON	2		
	MN	MANGANESE	TRACE		
	FE	IRON	TRACE		
	C	CARBON	TRACE		
Titanium	TI	TITANIUM	100	SILVER	IMPLANT ABUTMENTS PORCELAIN BONDING
Argenco 60	AU	GOLD	60	RICH YELLOW	CAST RESTORATIONS 60% precious metal Suitable for NHS Inlay/Onlay
	AG	SILVER	21.75		
	PD	PALLADIUM	3.5		
	CU	COPPER	13.49		
	ZN	ZINC	1		
	IN	INDIUM	TRACE		
	IR	IRIDIUM	TRACE		
Argenco 10S	AU	GOLD	10	PALE YELLOW	CAST RESTORATIONS 34% precious metal Suitable for NHS 3/4 Crowns
	PD	PALLADIUM	24		
	AG	SILVER	45.9		
	IN	INDIUM	18		
	ZN	ZINC	2		
	IR	IRIDIUM	TRACE		
NPG/Argeloy Sunray	CU	COPPER	80.7	RICH YELLOW	CAST RESTORATIONS non-precious metal Not suitable for NHS restorations
	AL	ALUMINUM	7.8		
	FE	IRON	3		
	ZN	ZINC	2.7		
	MN	MANGANESE	1.7		
	NI	NICKEL	4.1		
Verasoft	NI	NICKEL	53.6	SILVER	CAST RESTORATIONS non-precious metal Little used due to Nickel Content
	MN	MANGANESE	19.5		
	CR	CHROMIUM	14.5		
	CU	COPPER	9.5		
	AL	ALUMINUM	1.6		
	SI	SILICON	1.6		
	FE	IRON	TRACE		